


Summary

Name	GRAP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75791
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEAVAKFDFT ASGEDELSFH TGDVLKILSN QEEWFKAELG SQEGYVPKNF IDIQFPKWFH EGLSRHQAEN LLMGKEVGFF IIRASQSSPG DFSISVRHED DVQHFKVMRD NKGNYFLWTE KFPSLNKLVD YYRTNSISRQ KQIFLRDRTR EDQGHRGNSL DRRSQGGPHL SGAVGEEIRP SMNRKLSDHP PTLPLQQHQH QPQPPQYAPA PQLLQPPQQ RYLQHHHFHQ ERRGGSLDIN DGHC GTGLGS EMNAALMHRH HTDPVQLQAA GRVRWARALY DFEALDEL GFHSGEVVEV LDSSNPSWWT GRLHNKLGFL PANYVAPMTR

image



Background

Alternative Names	GRB2-related adapter protein 2 (Adapter protein GRID) (GRB-2-like protein) (GRB2L) (GRBLG) (GRBX) (Grf40 adapter protein) (Grf-40) (Growth factor
--------------------------	---

Product Name: Recombinant Human GRAP2 Protein
Catalog #: PEH4317



receptor-binding protein) (Hematopoietic cell-associated adapter protein
GrpL) (P38) (Protein GADS) (SH3-SH2-SH3 adapter Mona)

Note

For Research Use Only , Not for Diagnostic Use