


Product Name: Recombinant Human SFTP A2 Protein
Catalog #: PEH4331



Summary

Name	SFPA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWL1
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVKDVCVGSP GIPGTPGSHG LPGRDGRDGV KGDPGPPGPM GPPGETPCPP GNNGLPGAPG VPGERGEKGE AGERGPPGLP AHLDEELQAT LHDFRHQILQ TRGALS LQGS IMTVGEKVFS SNGQSITFDA IQEACARAGG RIAVPRNPEE NEIASFVKK YNTYAYVGLT EGSPGDFRY SDGTPVNYTN WYRGEPAGRG KEQCVEMYTD GQWNDNRNCLY SRLTICEF

image



Background

Alternative Names	Pulmonary surfactant-associated protein A2 (PSP-A) (PSPA) (SP-A) (SP-A2) (35 kDa pulmonary surfactant-associated protein) (Alveolar proteinosis protein) (Collectin-5)
--------------------------	--

Product Name: Recombinant Human SFTPA2 Protein
Catalog #: PEH4331



Note

For Research Use Only , Not for Diagnostic Use