


Product Name: Recombinant Human PTPN6 Protein
Catalog #: PEH4336



Summary

Name	PTN6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P29350
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPYATRVNA ADIENRVLEL NKKQESDTA KAGFWEFES LKQQEVKNLH QRLEGQRPEN KGKNRYKNIL PFDHSRVILQ GRDSNIPGSD YINANYIKNQ LLGPDENAKT YIASQGCLEA TVNDFWQMAW QENSRVIVMT TREVEKGRNK CVPYWPEVGM QRAYGPYSVT NCGEHDTEY KLRTLQVSPL DNGDLIREIW HYQYLSWPDH GVPSEPGGV L SFLDQINQRQ ESLPHAGPII VHCSAGIGRT GTIIVIDMLM ENISTKGLDC DIDIQKTIQM VRAQRSGMVQ TEAQYKFIYV AIAQFIETTK KKLEVLSQSK GQSEYGNIT YPPAMKNAHA KASRTSSKHK EDVYENLHTK NKREEKVKKQ RSADKEKSKG

image



Background

Alternative Names Tyrosine-protein phosphatase non-receptor type 6 (EC 3.1.3.48)

Product Name: Recombinant Human PTPN6 Protein
Catalog #: PEH4336



(Hematopoietic cell protein-tyrosine phosphatase) (Protein-tyrosine phosphatase 1C) (PTP-1C) (Protein-tyrosine phosphatase SHP-1) (SH-PTP1)

Note

For Research Use Only , Not for Diagnostic Use