



Product Name: Recombinant Human SEMA3F Protein
Catalog #: PEH4337



Summary

Name	SEM3F_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q13275
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVHAALFPPL SMSAPPPGA GPPTPPYQEL AQLLAQPEVG LIHQYCQGYW RHHVPPSPREA PGAPRSPEPQ

image



Background

Alternative Names Semaphorin-3F (Sema III/F) (Semaphorin IV) (Sema IV)

Note

For Research Use Only , Not for Diagnostic Use