

Product Name: Recombinant Human CES2 Protein
Catalog #: PEH4343



Summary

Name	EST2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00748
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLGSLVHVKG ANAGVQTF LG IPFAKPPLGP LRFAPPEPPE SWSGVRDGT T HPAMCLQDLT AVESEFLSQF NMTFPSDSMS EDCLYLSIYT PAHSHEGSNL PVMVWIHGGA LVFGMASLYD GSMLAALENV VVVIQYRLG VLGFFSTGDK HATGNWGYLD QVAALRWVQQ NIAHFGGNPD RVTIFGESAG GTSVSSLVVS PISQGLFHGA IMESGVALLP GLIASSADVI STVVANLSAC DQVDSEALVG CLRGKSKEEI LAINKPFKMI PGVVDGVFLP RHPQELLASA DFQPVPSIVG VNNNEFGWLI PKVMRIYDTQ KEMDREASQA ALQKMLTLLM LPPTFGDLLR EEYIGDNGDP QTLQAQFQEM MADSMFVIPA LQVAHFQCSR APVYFYEFQH QPSWLKNIRP PHMKADHGDE LPFVFRSFFG GNYIKFTEEE EQLSRKMMKY WANFARNGNP NGEGLPHWPL FDQEEQYLQL NLQPAVGRAL KAHRLQFWKK ALPQKIQELE

image



Background

Product Name: Recombinant Human CES2 Protein
Catalog #: PEH4343



Alternative Names

Cocaine esterase (EC 3.1.1.84) (Carboxylesterase 2) (CE-2) (hCE-2) (EC 3.1.1.1) (Methylumbelliferyl-acetate deacetylase 2) (EC 3.1.1.56)

Note

For Research Use Only , Not for Diagnostic Use