


Product Name: Recombinant Human ITGB1BP1 Protein
Catalog #: PEH4344



Summary

Name	ITBP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14713
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TKSKSV DSSL GGLSR SSTVA SLDTD STKSS GQSNN NSDTC AEFRI KYVGA IEK LK LSE GK GLEGL DLIN YIDVA QQDGK LPFVP PEEEF IMGVSK YGIK VSTSD QYDVL HRHAL YLIIR MVCYD DGLGA GKSL LALKTT DASNEE YSLW VYQCNS LEQA QAICK VLSTA FDSVLT SEKP

image



Background

Alternative Names	Integrin beta-1-binding protein 1 (Integrin cytoplasmic domain-associated protein 1) (ICAP-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ITGB1BP1 Protein
Catalog #: PEH4344

