


Product Name: Recombinant Human EPHB6 Protein
Catalog #: PEH4346



Summary

Name	EPHB6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O15197
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAPQELWFEV QGSALMLHWR LPRELGGRGD LLFNVVCKEC EGRQEPASGG GGTCHRCRDE VHFDPQRQRL TESRVLVGGL RAHVPYILEV QAVNGVSELS PDPPQAAAIN VSTSHEVPSA VPVVHQVSRA SNSITVSWPQ PDQTNGNILD YQLRYYDQAE DESHSFTLTS ETNTATVTQL SPGHIYGFQV RARTAAGHGP YGGKVFYFQTL

image



Background

Alternative Names	Ephrin type-B receptor 6 (HEP) (Tyrosine-protein kinase-defective receptor EPH-6)
--------------------------	---

Note

Product Name: Recombinant Human EPHB6 Protein
Catalog #: PEH4346



For Research Use Only , Not for Diagnostic Use