


Product Name: Recombinant Human FGF10 Protein
Catalog #: PEH4348



Summary

Name	FGF10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O15520
Host	E.coli
Species	Human
Molecular Mass	62 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSSVPVTCQA LGQDMVSPEA TNSSSSFSS PSSAGRHVRS YNHLQGDVRW RKLFSTKYF LKIEKNGKVS GTKKENCYPYS ILEITSVEIG VVAVKAINSN YYLAMNKKGK LYGSKEFNND CKLKERIEEN GYNTYASFNW QHNGRQMYVA LNGKGAPRRG QKTRRKNTSA HFLPMVVHS

image



Background

Alternative Names Fibroblast growth factor 10 (FGF-10) (Keratinocyte growth factor 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FGF10 Protein
Catalog #: PEH4348

