


Product Name: Recombinant Human DNPH1 Protein
Catalog #: PEH4351



Summary

Name	DNPH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43598
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAAMVPGRS ESWERGEPGR PALYFCGSIR GGREDRTLYE RIVSRLRRFG TVLTEHVAAA ELGARGEAAA GGDRLIHEQD LEWLQQADV VAEVTQPSLG VGYELGRAVA FNKRILCLFR PQSGRVLSAM IRGAADGSRF QVWDYEEGEV EALLDRYFEA DPPGQVAASP DPTT

image



Background

Alternative Names 2'-deoxynucleoside 5'-phosphate N-hydrolase 1 (EC 3.2.2.-) (c-Myc-responsive protein RCL)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DNPH1 Protein
Catalog #: PEH4351

