


Product Name: Recombinant Human XPNPEP2 Protein
Catalog #: PEH4354



Summary

Name	XPP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43895
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QIEDYSQVRD SIQAYSLGDV RIWIGTSYTM YGIYEMIPKE KLVTDITYSPV MMTKAVKNSK EQALLKASHV RDAVAVIRYL VWLEKNVPKG TVDEFSGAEI VDKFRGEEQF SSGPSFETIS ASGLNAALAH YSPTKELNRK LSSDEMYLLD SGGQYWDGTT DITRTVHWGT PSAFQKEAYT RVLIGNIDLS RLIFPAATSG RMVEAFARRA LWDAGLNYGH GTGHGIGNFL CVHEWPVGFQ SNNIAMAKGM FTSIEPGYYK DGEFGIRLED VALVVEAKTK

image



Background

Alternative Names	Xaa-Pro aminopeptidase 2 (EC 3.4.11.9) (Aminoacylproline aminopeptidase) (Membrane-bound aminopeptidase P) (Membrane-bound APP) (Membrane-bound Amp) (mAmP) (X-Pro aminopeptidase 2)
--------------------------	--

Product Name: Recombinant Human XPNPEP2 Protein
Catalog #: PEH4354



Note

For Research Use Only , Not for Diagnostic Use