


Product Name: Recombinant Human ENTPD6 Protein
Catalog #: PEH4358



Summary

Name	ENTP6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O75354
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LVSPCLSPSF KGEWEHAEVT YRVSGQKAAA SLHELCAARV SEVLQNRVHR TEEVKHFVDFY AFSYYDLAA GVGLIDAEKG GSLVVGDFEI AAKYVCRTLE TQPQSSPFSC MDLTYVLLLL QEFGFPRSKV LKLTRKIDNV ETSWALGAIF HYIDSLNRQK

image



Background

Alternative Names Ectonucleoside triphosphate diphosphohydrolase 6 (NTPDase 6) (EC 3.6.1.6) (CD39 antigen-like 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ENTPD6 Protein
Catalog #: PEH4358

