


Product Name: Recombinant Human CNTN5 Protein
Catalog #: PEH4362



Summary

Name	CNTN5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O94779
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGRSGRRHEL VIAWEPVSEE FQNGEGFGYI VAFRPNNGTRG WKEKMVTSSE ASKFIYRDES VPPLTPFEVK VGVYNNKGDG PFSQIVVICS AEGEPSAAPT DVKATSVSVS EILVAWKHIK ESLGRPQGFE VGYWKDMEQE DTAETVKTRG NESFVILTGL EGNTLYHFTV RAYNGAGYGP PSSEVSATTK KSPPSQAPSN LRWEQQGSQV SLGWEPVIPL ANESEVVGKY VFYRQEGHSN SQVIETQKLQ AVVPLPDAGV YIIEVRAYSE GGDGTASSQI

image



Background

Alternative Names	Contactin-5 (Neural recognition molecule NB-2) (hNB-2)
--------------------------	--

Product Name: Recombinant Human CNTN5 Protein
Catalog #: PEH4362



Note

For Research Use Only , Not for Diagnostic Use