


Product Name: Recombinant Human PLPHP Protein
Catalog #: PEH4363



Summary

Name	PLPHP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94903
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MWRAGSMSAE LGVGCALRAV NERVQQAVAR RPRDLPAIQP RLVAVSKTKP ADMVIEAYGH GQRTFGENYV QELLEKASNP KILSLCPEIK WHFIGHLQKQ NVNKLMAVPN LFMLETVDVSV KLADKVNSSW QRKGSPELRLK VMVQINTSGE ESKHGLPPSE TIAIVEHINA KCPNLEFVGL MTIGSFGHDL SQGPNPDFQL LLSLREELCK KLNIPADQVE LSMGMSADFQ HAVEVGSTNV RIGSTIFGER DYSKKPTPDK CAADVKAPEL VAQEH

image



Background

Alternative Names	Pyridoxal phosphate homeostasis protein (PLP homeostasis protein) (Proline synthase co-transcribed bacterial homolog protein) (Pyridoxal phosphate-binding protein)
--------------------------	---

Product Name: Recombinant Human PLPBP Protein
Catalog #: PEH4363



Note

For Research Use Only , Not for Diagnostic Use