


Product Name: Recombinant Human CST3 Protein
Catalog #: PEH4371



Summary

| | |
|------------------------------|--|
| Name | CYTC_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-SUMO & C-His |
| Accession # | P01034 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 26 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SSPGKPPRLV GGPMDSVVEE EGVRRALDFA VGEYNKASND MYHSRALQVV RARKQIVAGV NYFLDVELGR TTCTKTQPNL DNCPFHDQPH LKRKAFCFSQ IYAVPWQGTM TSKSTCQDA |

image



Background

Alternative Names Cystatin-C (Cystatin-3) (Gamma-trace) (Neuroendocrine basic polypeptide) (Post-gamma-globulin)

Note

For Research Use Only , Not for Diagnostic Use