


Summary

Name	ITAV_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P06756
Host	E.coli
Species	Human
Molecular Mass	76 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PCRSQDIDAD GQGFCQGGFS IDFTKADRVL LGGPGSFYWQ GQLISDQVAE IVSKYDPNVY SIKYNNQLAT RTAQAIFFDS YLGYSVAVGD FNGDGIDDFV SGVPRAARTL GMVYIYDGKN MSSLYNFTGE QMAAYFGFSV AATDINGDDY ADVFIGAPLF MDRGSDGKLQ EVGQVSVSLQ RASGDFQTTK LNGFEVFARF GSAIAPLGD L DQDGFNDIAI AAPYGGEDKK GIVYIFNGRS TGLNAVPSQI LEGQWAARSM PPSFGYSMKG ATDIDKNGYP DLIVGAFGVD RAILYRARPV ITVNAGLEVY PSILNQDNKT CSLPGTALKV SCFNVRFCLK ADGKGVLPK LNFQVELLLD KLKQKGAIRR ALFLYSRSPS HSKNMTISRG GLMQCEELIA YLRDESEFRD KLTPIIFME YRLDYRTAAD TTGLQPILNQ FTPANISRQA HILLDCGEDN VCKPKLEVS DSDQKKIYIG DDNPLTLIVK AQNQGEGAYE AELIVSIPLQ ADFIGVVRNN EALARLSCAF K TENQTRQVV CDLGNP MKAG TQLLAGLRFS VHQQSEM DTS

image



Product Name: Recombinant Human ITGAV Protein
Catalog #: PEH4384



Background

Alternative Names

Integrin alpha-V (Vitronectin receptor) (Vitronectin receptor subunit alpha) (CD antigen CD51) [Cleaved into: Integrin alpha-V heavy chain; Integrin alpha-V light chain]

Note

For Research Use Only , Not for Diagnostic Use