


Summary

| | |
|------------------------------|---|
| Name | MMP10_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-MBP & C-His |
| Accession # | P09238 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 85 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | FSSFPGMPKW RKTHLYRIV NYTPDLPRDA VDSAIEKALK VWEEVTPLTF SRLYEGEADI MISFAVKEHG DFYSFDGPGH SLAHAYPPGP GLYGDHIFDD DEKWTEASG TNLFLVAAHE LGHSLGLFHS ANTEALMYPL YNSFTELAQF RLSQDDVNGI QSLYGPPAS TEEPLVPTKS VPSGSEMPAK CDPALSFDAI STLRGEYLFF KDRYFWRRSH WNPEPEFHLL SAFWPSLPSY LDAAYEVNSR DTVFIKGFNE FWAIRGNEVQ AGYPRGIHTL GFPPTIRKID AAVSDKEKKK TYFFAADKYW RFDENSQSME QGFPRLIADD FPGVEPKVDA VLQAFGFFYF FSGSSQFEFD PNARMVTHIL KSNSWLHC |

image



Background

Alternative Names Stromelysin-2 (SL-2) (EC 3.4.24.22) (Matrix metalloproteinase-10) (MMP-10)

Product Name: Recombinant Human MMP10 Protein
Catalog #: PEH4394



(Transin-2)

Note

For Research Use Only , Not for Diagnostic Use