


Product Name: Recombinant Human QDPR Protein
Catalog #: PEH4396



Summary

Name	DHPR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09417
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAAAAGEA RRVLVYGGRG ALGSRVQAF RARNWWVASV DVVENEESA SIIVKMTDSF TEQADQVTAE VGKLLGEEKV DAILCVAGGW AGGNAKSKSL FKNCDLMWKQ SIWTSTISSH LTKHLKEGG LLTLAGAKAA LDGTPGMIGY GMAKGAVHQL CQSLAGKNSG MPPGAAAIIV LPVTLDTPMN RKSMPPEAFS SWTPLEFLVE TFHDWITGKN RPSSGSLIQV VTTEGRTELT PAYF

image



Background

Alternative Names Dihydropteridine reductase (EC 1.5.1.34) (HDHPR) (Quinoid dihydropteridine reductase) (Short chain dehydrogenase/reductase family 33C member 1)

Note

Product Name: Recombinant Human QDPR Protein
Catalog #: PEH4396



For Research Use Only , Not for Diagnostic Use