


Product Name: Recombinant Human CTSB Protein
Catalog #: PEH4398



Summary

Name	CATH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09668
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESAIAIATGK MSLAEQQLV DCAQDFNNHG CQGGLPSQAF EYILYNKGIM GEDTYPYQ GK DGYCKFQPGK AIGFVKDVAN ITIYDEEAMV EAVALYNPVS

image



Background

Alternative Names Pro-cathepsin H [Cleaved into: Cathepsin H mini chain; Cathepsin H (EC 3.4.22.16); Cathepsin H heavy chain; Cathepsin H light chain]

Note

For Research Use Only , Not for Diagnostic Use