


Product Name: Recombinant Human TACSTD2 Protein
Catalog #: PEH4400



Summary

Name	TACD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P09758
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	THHILIDLRH RPTAGAFNHS DLDAELRRLF RERYRLHPKF VAAVHYEQPT IQIELRQNTS QKAAGDVDIG DAAYYFERDI KGESLFQGRG GLDLRVRGEP LQVERTLIYY LDEIPPKFSM

image



Background

Alternative Names	Tumor-associated calcium signal transducer 2 (Cell surface glycoprotein Trop-2) (Membrane component chromosome 1 surface marker 1) (Pancreatic carcinoma marker protein GA733-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TACSTD2 Protein
Catalog #: PEH4400

