


Product Name: Recombinant Human FGR Protein
Catalog #: PEH4401



Summary

Name	FGR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P09769
Host	E.coli
Species	Human
Molecular Mass	70 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVATAKEDAG LEGDFRSYGA ADHYGPDPTK ARPASSFAHI PNYSNFSSQA INPGFLDSGT IRGVSGIGVT LFIALYDYEY RTEDDLTFTK GEKFHILNNT EGDWWEARSL SSGKTGCIPS NYVAPVDSIQ AEEWYFGKIG RKDAERQLLS PGNPQGAFIL RESETTKGAY SLSIRDWDQT RGDHVKHYKI RKLDMGGYYI TTRVQFNSVQ ELVQHMEVN DGLCNLLIAP CTIMKPQTLG LAKDAWEISR

image



Background

Alternative Names	Tyrosine-protein kinase Fgr (EC 2.7.10.2) (Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog) (Proto-oncogene c-Fgr) (p55-Fgr) (p58-Fgr) (p58c-Fgr)
--------------------------	---

Product Name: Recombinant Human FGR Protein
Catalog #: PEH4401



Note

For Research Use Only , Not for Diagnostic Use