


Product Name: Recombinant Human FGR Protein
Catalog #: PEH4402



Summary

Name	FGR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09769
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RTEDDLTFTK GEKFHILNNT EGDWWEARSL SSGKTGCIPS NYVAPVDSIQ AEEWYFGKIG RKDAERQLLS PGNPQGAFLI RESETTKGAY SLSIRDWDQT RGDHVKHYKI RKLDMMGGYYI TTRVQFNSVQ ELVQHMEVN

image



Background

Alternative Names	Tyrosine-protein kinase Fgr (EC 2.7.10.2) (Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog) (Proto-oncogene c-Fgr) (p55-Fgr) (p58-Fgr) (p58c-Fgr)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FGR Protein
Catalog #: PEH4402

