



Product Name: Recombinant Human LTA4H Protein
Catalog #: PEH4404



Summary

Name	LKHA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P09960
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CRAILPCQDT PSVKLTYTAE VSVPKELVAL MSAIRDGETP DPEDPSRKIY KFIQKVPIPC YLIALVVGAL ESRQIGPRTL VWSEKEQVEK SAYEFSETES MLKIAEDLGG PYVWGQYDLL VLPPSFPYGG MENPCLTFVT PTLLAGDKSL SNVIAHEISH SWTGNLVTNK TWDHFWLNEG HTVYLERHIC GRLFGKEFRH FNALGGWGEL QNSVKTFGET HPFTKLVDL TDIDPDVAYS SVPYEKGFLAL LFYLEQLLGG PEIFLGFLKA

image



Background

Alternative Names	Leukotriene A-4 hydrolase (LTA-4 hydrolase) (EC 3.3.2.6) (Leukotriene A(4) hydrolase) (Tripeptide aminopeptidase LTA4H) (EC 3.4.11.4)
--------------------------	---

Product Name: Recombinant Human LTA4H Protein
Catalog #: PEH4404



Note

For Research Use Only , Not for Diagnostic Use