


Product Name: Recombinant Human CKMT1A Protein
Catalog #: PEH4409



Summary

Name	KCRU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12532
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEYPLRKHN NCMASHLTPA VYARLCDKTT PTGWTLDQCI QTGVDNPGHP FIKTVGMVAG DEETYEVFAD LFDPVIQERH NGYDPRTMKH TTDLDASKIR SGYFDERYVL SSRVRTGRSI RGLSLPPACT RAERREVERV VVDALSGLK DLAGRYRRLS EMTEAEQQQL IDHFLFDKP VSPLLTAAGM ARDWPDARGI WHNNEKSFLI WVNEEDHTRV ISMEKGGNMK RVFERFCRGL KEVERLIQER GWEFMWNERL GYILTCPSNL GTGLRAGVHI KLPLSKDSR FPKILENLRL QKRGTTGGVDT AATGGVFDIS NLDRLGKSEV ELVQLVIDGV NYLIDCERRL ERGQDIRIPT

image



Background

Alternative Names	Creatine kinase U-type, mitochondrial (EC 2.7.3.2) (Acidic-type mitochondrial
--------------------------	---

Product Name: Recombinant Human CKMT1A Protein
Catalog #: PEH4409



creatine kinase) (Mia-CK) (Ubiquitous mitochondrial creatine kinase) (U-MtCK)

Note

For Research Use Only , Not for Diagnostic Use