


Product Name: Recombinant Human PEPD Protein
Catalog #: PEH4412



Summary

Name	PEPD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P12955
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CERLRKNPAV QAGSIVVLQG GEETQRYCTD TGVLFQESF FHWAFGVTEP GCYGVIDVDT GKSTLFVPRP PASHATWMGK IHSKEHFKEK YAVDDVQYVD EIASVLTSQK

image



Background

Alternative Names	Xaa-Pro dipeptidase (X-Pro dipeptidase) (EC 3.4.13.9) (Imidodipeptidase) (Peptidase D) (Proline dipeptidase) (Prolidase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use