


Product Name: Recombinant Human ARSA Protein
Catalog #: PEH4417



Summary

Name	ARSA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15289
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSGRGPFQDS LMEIDAAVGT LMTAIGDLGL LEETLVIFTA DNGPETMRMS RGGCSGLLRC GKGTTEGGV REPALAFWPG HIAPGVTHEL ASSLDLLPTL AALAGAPLPN VTLDGFDLSP LLLGTGKSPR QSLFFYPSTP DEVRGVFAVR TGKYKAHFFT QGSAHSDDTA DPACHASSL TAHEPPLLYD LSKDPGENYN LLGGVAGATP EVLQALKQLQ LLKAQLDAV TFGPSQVARG

image



Background

Alternative Names	Arylsulfatase A (ASA) (EC 3.1.6.8) (Cerebroside-sulfatase) [Cleaved into: Arylsulfatase A component B; Arylsulfatase A component C]
--------------------------	---

Note

Product Name: Recombinant Human ARSA Protein
Catalog #: PEH4417



For Research Use Only , Not for Diagnostic Use