


Product Name: Recombinant Human ARSB Protein
Catalog #: PEH4418



Summary

| | |
|------------------------------|--|
| Name | ARSB_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P15848 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 28 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | RCALDFRDGE EVATGYKNMY STNIFTKRAI ALITNHPPEK PLFLYLALQS VHEPLQVPEE YLKPYDFIQD KNRHHYAGMV SLMDEAVGNV TAALKSSGLW NNTVFIFSTD NGGQTLAGGN NWPLRGRKWS LWEGGVRGVG FVASPLLKQK GVKNRELIHI SDWLPTLVKL ARGHTNGTKP LDGFDVWKTI SEGSPSPRIE LLHNIDPNFV DSSPCPRNSM APAKDDSSLP |

image



Background

Alternative Names Arylsulfatase B (ASB) (EC 3.1.6.12) (N-acetylgalactosamine-4-sulfatase) (G4S)

Note

Product Name: Recombinant Human ARSB Protein
Catalog #: PEH4418



For Research Use Only , Not for Diagnostic Use