


**Product Name: Recombinant Human ARSB Protein**  
**Catalog #: PEH4419**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | ARSB_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | N-SUMO & C-His   |
| <b>Accession #</b>           | P15848   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 39 kDa   |
| <b>Buffer</b>                | Lyophilized from Phosphate buffered saline   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.  |
| <b>Sequence</b>              | RCALDFRDGE EVATGYKNMY STNIFTKRAI ALITNHPPEK PLFLYLALQS<br>VHEPLQVPEE YLKPYDFIQD KNRHHYAGMV SLMDEAVGNV TAALKSSGLW<br>NNTVFIFSTD NGGQTLAGGN NWPLRGRKWS LWEGGVRGVG FVASPLLKQK<br>GVKNRELIHI SDWLPTLVKL ARGHTNGTKP LDGFDVWKTI SEGSPSPRIE<br>LLHNIDPNFV DSSPCPRNSM APAKDDSSLP |

## image



## Background

**Alternative Names** Arylsulfatase B (ASB) (EC 3.1.6.12) (N-acetylgalactosamine-4-sulfatase) (G4S)

## Note

**Product Name: Recombinant Human ARSB Protein**  
**Catalog #: PEH4419**



---

For Research Use Only , Not for Diagnostic Use