


Product Name: Recombinant Human COL9A1 Protein
Catalog #: PEH4422



Summary

Name	CO9A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20849
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FPVNSNSNGG NELCPKIRIG QDDLPGFDLI SQFQVDKAAS RRAIQRVVG ATLQVAYKLG NNVDFRIPTR NLYPSGLPEE YSFLTFRMT GSTLKKNWNI WQIQDSSGKE QVGIKINGQT QSVVFSYKGL DGSLQTAAFS NLSSLFDSQW HKIMIGVERS SATLFVDCNR IESLPIKPRG PIDIDGFAVL GKLADNPQVS VPFELQWMLI HCDPLRPRE TCHELPARIT PSQTTDERGP PGEQGP GPPGVPIDG IDGDRGPKGP PGPPGPAGEP GKPGAPGKPG

image



Background

Alternative Names Collagen alpha-1(IX) chain

Product Name: Recombinant Human COL9A1 Protein
Catalog #: PEH4422



Note

For Research Use Only , Not for Diagnostic Use