



## Summary

<b>Name</b>	RASA1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P20936
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	68 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	IPLTAPPTNQ WYHGKLDRTI AEERLRQAGK SGSYLIRESD RRPGSFVLSF LSQMNVVNHF RIIAMCGDYY IGGRRFSSLS DLIGYYSHVS CLLKGEKLLY PVAPPEPVED RRRVRAILPY TKVPDTDEIS FLKGDMFIVH NELEDGWMWV TNLRTDEQGL IVEDLVEEVG REEDPHEGKI WFHGKISKQE AYNLLMTVGQ VCSFLVRPSD NTPGDYSLYF RTNENIQRFK ICPTPNNQFM MGGRYNSIG DIIDHYRKEQ IVEGYLKEP VPMQDQEQVL NDTV DGKEIY NTIRRKTKDA FYKNIVKKG Y LLKKGKGRW KNLYFILEGS DAQLIFYESE KRATKPKGLI DLSVCSVYV HDSLFRPNC FQIVVQHSE EHYIFYFAGE TPEQAEDWMK GLQAF CNLRK SSPGTSNKRL RQVSSLV LHI EEAHKL PVKH FTNPYCN IYL NSVQVAKTHA REGQNPVWSE EFVFDLPPD INRFEITLSN KTKKSKDPDI LFMRCQLSRL QKGHATDEWF LLSSHIPLKG IEPGSLRVRA RYSMEKIMPE EEYSEFKELI LQKELHV VYA

## image



**Product Name: Recombinant Human RASA1 Protein**  
**Catalog #: PEH4423**



---

## Background

**Alternative Names**

Ras GTPase-activating protein 1 (GAP) (GTPase-activating protein) (RasGAP)  
(Ras p21 protein activator) (p120GAP)

## Note

For Research Use Only , Not for Diagnostic Use