



Product Name: Recombinant Human RASA1 Protein
Catalog #: PEH4424



Summary

Name	RASA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P20936
Host	E.coli
Species	Human
Molecular Mass	80 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IPLTAPPTNQ WYHGKLDRTI AEERLRQAGK SGSYLIRESD RRPGSFVLSF LSQMNVVNHF RIIAMCGDYY IGRRFSSLS DLIGYYSHVS CLLKGEKLLY PVAPPEPVED RRRVRAILPY TKVPDTDEIS FLKGDMFIVH NELEDGWMWV TNLRTDEQGL IVEDLVEEVG REEDPHEGKI WFHGKISKQE AYNLLMTVGQ VCSFLVRPSD NTPGDYSLYF RTNENIQRFK ICPTPNNQFM MGGRYNSIG DIIDHYRKEQ IVEGYLKEP VPMQDQEQVL NDTV DGKEIY NTIRRKT KDA FYKNIVKKG Y LLKKGKGKRW KNLYFILEGS DAQLIFYESE KRATKPKGLI DLSVCSVYV HDSLFGRPNC FQIVVQH FSE EHYIFYFAGE TPEQAEDWMK GLQAF CNLRK SSPGTSNKRL RQVSSLV LHI EEAHKLPVKH FTNPYCNIYL NSVQVAKTHA REGQNPVWSE EFVFDDLPPD INRFEITLSN KTKKSKDPDI LFMRCQLSRL QKGHATDEWF LLSSHIPLKG IEPGSLRVRA RYSMEKIMPE EEYSEFKELI LQKELHV VYA

image



Product Name: Recombinant Human RASA1 Protein
Catalog #: PEH4424



Background

Alternative Names

Ras GTPase-activating protein 1 (GAP) (GTPase-activating protein) (RasGAP)
(Ras p21 protein activator) (p120GAP)

Note

For Research Use Only , Not for Diagnostic Use