


Product Name: Recombinant Human NT5E Protein
Catalog #: PEH4425



Summary

Name	5NTD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P21589
Host	E.coli
Species	Human
Molecular Mass	73 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKYPFIVTSD DGRKVPVVQA YAFGKYLGYL KIEFDERGNV ISSHGNPILL NSSIPEDPSI KADINKWRIK LDNYSTQELG KTIVYLDGSS QSCRFRECNM GNLICDAMIN NNLRHTDEMF WNHVSMCILN GGGIRSPIDE RNNGTITWEN LAAVLPGGT FDLVQLKGST LKKA FEHSVH RYQSTGEFL QVGGIHVVYD LSRKPGDRVV KLDVLCTKCR VPSYDPLKMD EVYKVILPNF LANGGDGFQM IKDELLRHDS

image



Background

Alternative Names 5'-nucleotidase (5'-NT) (EC 3.1.3.5) (Ecto-5'-nucleotidase) (CD antigen CD73)

Product Name: Recombinant Human NT5E Protein
Catalog #: PEH4425



Note

For Research Use Only , Not for Diagnostic Use