


Product Name: Recombinant Human NT5E Protein
Catalog #: PEH4426



Summary

Name	5NTD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P21589
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKYPFIVTSD DGRKVPVVQA YAFGKYLGYL KIEFDERGNV ISSHGNPILL NSSIPEDPSI KADINKWRIK LDNYSTQELG KTIVYLDGSS QSCRFRECNM GNLICDAMIN>NNLRHTDEMF WNHVSMCILN GGGIRSPIDE RNNGTITWEN LAAVLPGGT FDLVQLKGST LKKA FEHSVH RYQSTGEFL QVGGIHVVYD LSRKPGDRVV KLDVLCTKCR VPSYDPLKMD EVYKVILPNF LANGGDGFQM IKDELLRHDS

image



Background

Alternative Names 5'-nucleotidase (5'-NT) (EC 3.1.3.5) (Ecto-5'-nucleotidase) (CD antigen CD73)

Product Name: Recombinant Human NT5E Protein
Catalog #: PEH4426



Note

For Research Use Only , Not for Diagnostic Use