


Product Name: Recombinant Human BGN Protein
Catalog #: PEH4429



Summary

Name	PGS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P21810
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASGADTSGVL DPDSVTPTYs AMCPFGCHCH LRVVQCSDLG LKSVPEISP DTLLDLQNN DISELRKDDF KGLQHLYALV LVNNKISKIH EKAFSPRLKL QKLYISK NHL VEIPP NLPSS LVELRIHDNR IRKVPKGVFS GLRNMNCIEM GGNPLENSGF EPGA FDGLKL NYLRISEAKL TGIPKDLPET LNELHLDHNK IQAI ELEDLL RYSKLYRLGL GHNQIRMIEN GSLSFLPTLR ELHLDNNKLA

image



Background

Alternative Names Biglycan (Bone/cartilage proteoglycan I) (PG-S1)

Note

Product Name: Recombinant Human BGN Protein
Catalog #: PEH4429



For Research Use Only , Not for Diagnostic Use