


Product Name: Recombinant Human GPX3 Protein
Catalog #: PEH4435



Summary

Name	GPX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P22352
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENSEILPTLK YVRPGGGFVP NFQLFEKGDV NGEKEQKFYT FLKNSCPPTS ELLGTSDFLF WEPMKVHDIR WNFEEKFLVGP DGIPIMRWHH RTTVSNVKMD ILSYMRRQAA

image



Background

Alternative Names Glutathione peroxidase 3 (GPx-3) (GSHPx-3) (EC 1.11.1.9) (Extracellular glutathione peroxidase) (Plasma glutathione peroxidase) (GPx-P) (GSHPx-P)

Note

For Research Use Only , Not for Diagnostic Use