



Product Name: Recombinant Human CES1 Protein
Catalog #: PEH4436



Summary

Name	EST1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P23141
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAPTMYEFQ YRPSFSDMK PKTVIGDHGD ELFSVFGAPF LKEGASEEEI RLSKMVMKFW ANFARNGNPN GEGLPHWPEY NQKEGYLQIG ANTQAAQKLK DKEVAFWTNL FAKKAVEKPP

image



Background

Alternative Names	Liver carboxylesterase 1 (Acyl-coenzyme A:cholesterol acyltransferase) (ACAT) (Brain carboxylesterase hBr1) (Carboxylesterase 1) (CE-1) (hCE-1) (EC 3.1.1.1) (Cholesteryl ester hydrolase) (CEH) (EC 3.1.1.13) (Cocaine carboxylesterase) (Egasyn) (HMSE) (Methylumbelliferyl-acetate deacetylase 1) (EC 3.1.1.56) (Monocyte/macrophage serine esterase) (Retinyl ester hydrolase) (REH) (Serine esterase 1) (Triacylglycerol hydrolase) (TGH)
--------------------------	--

Product Name: Recombinant Human CES1 Protein
Catalog #: PEH4436



Note

For Research Use Only , Not for Diagnostic Use