


Product Name: Recombinant Human IL32 Protein
Catalog #: PEH4443



Summary

Name	IL32_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P24001
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVGHLPWRD KDPALWCQLC LSSQHQAIEF FYDKMQNAES GRGQVMSSLA ELEDDFKEGY LETVAAYYEE QHPELTPLLE KERDGLRCRG NRSPVPDVED PATEEPGESF CDKVMRWFOA MLQRLQTWWH GVLAWVKEKV

image



Background

Alternative Names Interleukin-32 (IL-32) (Natural killer cells protein 4) (Tumor necrosis factor alpha-inducing factor)

Note

For Research Use Only , Not for Diagnostic Use