


Product Name: Recombinant Human CTSS Protein
Catalog #: PEH4445



Summary

Name	CATS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P25774
Host	E.coli
Species	Human
Molecular Mass	66 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPDSVDWREK GCVTEVKYQG SCGACWAFSA VGALEAQLKL KTGKLVLSA QNLVDCSTEK YGNKGCNGGF MTTAFQYIID NKGIDSDASY PYKAMDQKCC YDSKYRAATC SKYTELPYGR EDVLKEAVAN KGPVSVGVD A RHPSFFLYRS GVYYEPSCTQ NVNHGVLVVG YGDLNGKEYW LVKNSWGHNF GEEGYIRMAR NKGNHCGIAS FPSYPEI

image



Background

Alternative Names Cathepsin S (EC 3.4.22.27)

Note

Product Name: Recombinant Human CTSS Protein
Catalog #: PEH4445



For Research Use Only , Not for Diagnostic Use