


**Product Name: Recombinant Human ITGB7 Protein**  
**Catalog #: PEH4447**



## Summary

<b>Name</b>	ITB7_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P26010
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	41 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ASGEAEARRC ARREELLARG CPLEELEPR GQQEVLDQDP LSQGARGEGA TQLAPQRVRV TLRPGEPQQL QVRFLRAEGY PVDLYYIMDL SYSMKDDLER VRQLGHALLV RLQEVTHSVR IGFGSFVDKT VLPFVSTVPS KLRHPCPTL ERCQSPFSFH HVLSLTGDAQ AFEREVGRQS VSGNLDSPG GFDAILQAAL CQEQIGWRNV SLLVFTSDD TFHTAGDGKL GGIFMPDGH CHLDSNGLYS RSTFDYPSV GQVAQALSAA NIQPIFAVTS AALPVYQELS KLIPKSAVGE LSEDSSNVVQ LIMDAYNSLS STVTLEHSSL PPGVHISYES QCEGPEKREG

## image



## Background

**Alternative Names** Integrin beta-7 (Gut homing receptor beta subunit)

**Product Name: Recombinant Human ITGB7 Protein**  
**Catalog #: PEH4447**

---



**Note**

For Research Use Only , Not for Diagnostic Use