


Product Name: Recombinant Human ITGB7 Protein
Catalog #: PEH4448



Summary

Name	ITB7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P26010
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASGEAEARRC ARREELLARG CPLEELEPR GQQEVLDQDP LSQGARGEGA TQLAPQRRVRV TLRPGEPQQL QVRFLRAEGY PVDLYYIMDL SYSMKDDLER VRQLGHALLV RLQEVTHSVR IGFGSFVDKT VLPFVSTVPS KLRHPCPTL ERCQSPFSFH HVLSLTGDAQ AFEREVGRQS VSGNLDSPG GFDAILQAAL CQEQIGWRNV SLLVFTSDD TFHTAGDGKL GGIFMPDGH CHLDSNGLYS RSTFDYPSV GQVAQALSAA NIQPIFAVTS AALPVYQELS KLIPKSAVGE LSEDSSNVVQ LIMDAYNSLS STVTLEHSSL PPGVHISYES QCEGPEKREG

image



Background

Alternative Names Integrin beta-7 (Gut homing receptor beta subunit)

Product Name: Recombinant Human ITGB7 Protein
Catalog #: PEH4448



Note

For Research Use Only , Not for Diagnostic Use