


Product Name: Recombinant Human IVD Protein
Catalog #: PEH4449



Summary

Name	IVD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P26440
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AINGLSEEQR QLRQTMAKFL QEHLAPKAQE IDRSNEFKNL REFWKQLGNL GVLGITAPVQ YGGSGGLGYL HVLVMEEISR ASGAVGLSYG AHSNLCINQL VRNGNEAQKE KYLPKLISGE YIGALAMSEP NAGSDVVSMK LKAEKKGNHY ILNGNKFWIT NGPDADVLIV YAKTDLAAPV ASRGITAFIV EKGMPGFSTS KKLDKLGMRG SNTCELIFED CKIPAANILG HENKGVYVLM SGLDLERLVL AGGPLGLMQA VLDHTIPYLH VREAFGQKIG HFQLMQGKMA DMYTRLMACR QYVYNVAKAC DEGHCTAKDC AGVILYSAEC ATQVALDGIQ CFGGNGYIND FPMGRFLRDA KLYEIGAGTS EVRRLVIGRA

image



Background

Alternative Names	Isovaleryl-CoA dehydrogenase, mitochondrial (IVD) (EC 1.3.8.4) (Butyryl-CoA
--------------------------	---

Product Name: Recombinant Human IVD Protein
Catalog #: PEH4449



dehydrogenase) (EC 1.3.8.1)

Note

For Research Use Only , Not for Diagnostic Use