


Product Name: Recombinant Human CNTF Protein
Catalog #: PEH4450



Summary

Name	CNTF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P26441
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IWLARKIRSD LTALTESYVK HQGLNKNINL DSADGMPVAS TDQWSELTEA ERLQENLQAY RTFHVLLARL LEDQQVHFTP TEGDFHQAIH TLLQVAAFA YQIEELMILL EYKIPRNEAD GMPINVGDDG LFEKKLWGLK VLQELSQWTV RSIHDLRFIS SHQTGIPARG

image



Background

Alternative Names Ciliary neurotrophic factor (CNTF)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CNTF Protein
Catalog #: PEH4450

