


**Product Name: Recombinant Human PON1 Protein**  
**Catalog #: PEH4452**



## Summary

<b>Name</b>	PON1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P27169
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LREVQPVELP NCNLVKGIET GSEDLEILPN GLAFISSGLK YPGIKSFNPN SPGKILLMDL NEEDPTVLEL GITGSKFDVS SFNPHGISTF TDEDNAMYLL VVNHPDAKST VELFKFQEEE KSLHLKTIR HKLLPNLNDI VAVGPEHFYQ TNDHYFLDPY LQSWEMYLGL AWSYVVYYSP SEVRVVAEGF DFANGINISP DGKYVYIAEL LAHKIHVYEK HANWTLPLK SLDFNTLVND ISVDPETGDL WVGCHPNGMK IFFYDSENP ASEVLRIQNI LTEEPKVTQV YAENGTVLQG STVASVYKGG LLIGTVFHKA

## image



## Background

**Alternative Names** Serum paraoxonase/arylesterase 1 (PON 1) (EC 3.1.1.2) (EC 3.1.1.81) (EC 3.1.8.1) (Aromatic esterase 1) (A-esterase 1) (K-45) (Serum

**Product Name: Recombinant Human PON1 Protein**  
**Catalog #: PEH4452**



---

aryldialkylphosphatase 1)

**Note**

For Research Use Only , Not for Diagnostic Use