



Product Name: Recombinant Human HSD11B1 Protein
Catalog #: PEH4457



Summary

Name	DHI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P28845
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AYHLAKMGAH VVVTARSKET LQKVVSHCLE LGAASAHYIA GTMEDMTFAE QFVAQAGKLM GGLDMLILNH ITNTSLNLFH DDIHHVRKSM EVNFLSYVVL TVAALPMLKQ SNGSIVVVSS LAGKVAYPMV AAYSASKFAL DGFFSSIRKE YSVSRVNVSI TLCVLGLIDT ETAMKAVSGI VHMQAAPKEE CALEIIKGGA

image



Background

Alternative Names	Corticosteroid 11-beta-dehydrogenase isozyme 1 (EC 1.1.1.146) (11-beta-hydroxysteroid dehydrogenase 1) (11-DH) (11-beta-HSD1) (Short chain dehydrogenase/reductase family 26C member 1)
--------------------------	---

Note

Product Name: Recombinant Human HSD11B1 Protein
Catalog #: PEH4457



For Research Use Only , Not for Diagnostic Use