


Product Name: Recombinant Human EPHA2 Protein
Catalog #: PEH4460



Summary

Name	EPHA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P29317
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLDFAAAGGE LGWLTHPYGK GWDLMQNIMN DMPIYMYSVC NVMSGDQDNW LRTNWVYRGE AERIFIELKF TVRDCNSFPG GASSCKETFN LYAESDL DY GTNFQKRLFT KIDTIAPDEI TVSSDFEARH VKLNVEERSV GPLTRKGFYL AFQDIGACVA LLSVRVYYKK CPELLQGLAH FPETIAGSDA PSLATVAGTC VDHAVVPPGG EEPRMHCAVD

image



Background

Alternative Names	Ephrin type-A receptor 2 (EC 2.7.10.1) (Epithelial cell kinase) (Tyrosine-protein kinase receptor ECK)
--------------------------	--

Note

Product Name: Recombinant Human EPHA2 Protein
Catalog #: PEH4460



For Research Use Only , Not for Diagnostic Use