


Product Name: Recombinant Human PRDX3 Protein
Catalog #: PEH4464



Summary

Name	PRDX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30048
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APAVTQHAPY FKGTAVVNGE FKDSLDDFK GKYLVLFFYP LDFTFVCPTE IVAFSDKANE FHDVNCEVVA VSVDSHFSL AWINTPRKNG GLGHMNIALL SDLTKQISRD YGVLLGSGSL ALRGLFIIDP NGVIKHLNVN DLPVGRSVEE TLRLVKAQFY VETHGEVCPA NWTPDSPTIK PSPAASKEYF QKVNQ

image



Background

Alternative Names	Thioredoxin-dependent peroxide reductase, mitochondrial (EC 1.11.1.24) (Antioxidant protein 1) (AOP-1) (HBC189) (Peroxiredoxin III) (Prx-III) (Peroxiredoxin-3) (Protein MER5 homolog) (Thioredoxin-dependent peroxiredoxin 3)
--------------------------	--

Product Name: Recombinant Human PRDX3 Protein
Catalog #: PEH4464



Note

For Research Use Only , Not for Diagnostic Use