


Product Name: Recombinant Human S100A11 Protein
Catalog #: PEH4470



Summary

Name	S10AB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-mSUMO & C-His
Accession #	P31949
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERCIESLIAV FQKYAGKDG Y NYTL SKTEFL SFMNTELA AF TKNQKDPGVL DRMMKKLDTN SDGQLDFSEF LNLIGGLAMA

image



Background

Alternative Names	Protein S100-A11 (Calgizzarin) (Metastatic lymph node gene 70 protein) (MLN 70) (Protein S100-C) (S100 calcium-binding protein A11) [Cleaved into: Protein S100-A11, N-terminally processed]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use