


Summary

Name	TAGLN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P37802
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MANRGPAYGL SREVQQKIEK QYDADLEQIL IQWITTQCRK DVGRPQPGRE NFQNWPKDGT VLCELINALY PEGQAPVKKI QASTMAFKQM EQISQFLQAA ERYGINTTDI FQTVDLWEGK NMACVQRTLM NLGGLAVARD DGLFSGDPNW FPKKSLENPR NFSDNQLQEG KNVIGLQMGT NRGASQAGMT GYGMPRQIL

image



Background

Alternative Names Transgelin-2 (Epididymis tissue protein Li 7e) (SM22-alpha homolog)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TAGLN2 Protein
Catalog #: PEH4475

