


Product Name: Recombinant Human CASP2 Protein
Catalog #: PEH4478



Summary

Name	CASP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P42575
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AYRLQSRPRG LALVLSNVHF TGEKELEFRS GGDVDHSTLV TLFKLLGYDV HVLCDQTAQE MQEKLQNFAQ LPAHRVTDSC IVALLSHGVE GAIYGVDGKL LQLQEVFQLF DNANCPSLQN KPKMFFIQAC

image



Background

Alternative Names	Caspase-2 (CASP-2) (EC 3.4.22.55) (Neural precursor cell expressed developmentally down-regulated protein 2) (NEDD-2) (Protease ICH-1) [Cleaved into: Caspase-2 subunit p18; Caspase-2 subunit p13; Caspase-2 subunit p12]
--------------------------	--

Note

Product Name: Recombinant Human CASP2 Protein
Catalog #: PEH4478



For Research Use Only , Not for Diagnostic Use