


Product Name: Recombinant Human PTPN9 Protein
Catalog #: PEH4479



Summary

Name	PTN9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P43378
Host	E.coli
Species	Human
Molecular Mass	75 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGIYEEYEDI RRENPGVTFH CSMSPGNLEK NRYGDVPCLD QTRVKLTKRS GHTQTDYINA SFMDGYKQKN AYIGTQGPLe NTYRDFWLMV WEQKVLVIVM TTRFEEGGRR KCGQYWPLEK DSRIRFGFLT VTNLGVENMN HYKKTLEIH NTEERQKRQV THFQFLSWPD YGVPSSAASL IDFLRVVRNQ QSLAVSNMGA RSKGQCPEPP IVVHCSAGIG RTGTFCSLDI CLAQLEELGT LNVFQTVSRM RTQRAFSIQT PEQYYFCYKA

image



Background

Alternative Names	Tyrosine-protein phosphatase non-receptor type 9 (EC 3.1.3.48) (Protein-tyrosine phosphatase MEG2) (PTPase MEG2)
--------------------------	--

Product Name: Recombinant Human PTPN9 Protein
Catalog #: PEH4479



Note

For Research Use Only , Not for Diagnostic Use