


Product Name: Recombinant Human SULT1A1 Protein
Catalog #: PEH4484



Summary

Name	ST1A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P50225
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLEYVKGVPL IKYFAEALGP LQSFQARPDD LLISTYPKSG TTWVSQILD IYQGGDLEKC HRAPIFMRVP FLEFKAPGIP SGMETLKDTP APRLLKTHLP LALLPQTLLD QKVKKVYVAR NAKDVAVSYY HFYHMAKVHP EPGTWDSFLE KFMVGEVSYG SWYQHVQEWV ELSRTHPVLY LFYEDMKENP KREIQKILEF VGRSLPEETV DFVVQHTSFK EMKKNPMTNY TTVPQEFMDH SISPFMRKGM AGDWKTFTV AQNERFDADY AEKMAGCSLS

image



Background

Alternative Names	Sulfotransferase 1A1 (ST1A1) (EC 2.8.2.1) (Aryl sulfotransferase 1) (HAST1/HAST2) (Phenol sulfotransferase 1) (Phenol-sulfating phenol sulfotransferase 1) (P-PST 1) (ST1A3) (Thermostable phenol sulfotransferase)
--------------------------	---

Product Name: Recombinant Human SULT1A1 Protein
Catalog #: PEH4484



(Ts-PST)

Note

For Research Use Only , Not for Diagnostic Use